


FDM Raman Retail Adhesives & Sealants



Samples:

3M Auto Plastic Emblem and Trim Adhesive
3M Carpenters Wood Glue, Indoor
3M Ceramic Adhesive, Outdoor, cured...

ELMERS(r) China plus Glass Cement
ELMERS(r) Contact Cement
ELMERS(r) Glue-All(r) Max...

Gold Eagle(r) AlumAseal(r) Radiator Stop Leak
Goo Gone(r) Pro Power(tm), PNGG12DO
Gorilla Glue, Precision Glue Pen

Features:

- 780 nm laser, 3250 cm⁻¹ to 200 cm⁻¹
- White Light Intensity Correction.
- Complements the FDM Raman Polymers & FDM ATR Retail Adhesives & Sealants.
- Spectral resolution check consistent
- with ASTM E2529-06.
- X axis calibration check consistent with ASTM E1840-96 (2007).
- 350 spectra run with new samples in 2013



Fiveash Data Management, Inc.

608-236-9145

www.fdmspectra.com

Fax: 608-236-9170