

Discover FDM's Very Large Bundle

Q: What is the FDM Very Large Bundle?

A: It's the best value in FTIR libraries.


Q: What formats are supported?

A: Agilent, Bruker, Essential FTIR, Jasco, PerkinElmer, Shimadzu, Thermo, Varian and others.

Q: How many libraries are included?

A: These 21 libraries, totaling more than 29,000 spectra.

- FDM ATR Polymers Diamond
- FDM ATR Polymers Germanium
- FDM ATR Organics Diamond
- FDM ATR Organics Germanium
- FDM ATR Inorganics Diamond
- FDM ATR Inorganics Germanium
- FDM ATR Dyes Diamond
- FDM ATR Dyes Germanium
- FDM ATR Cosmetics
- FDM ATR Retail Adhesives & Sealants
- FDM ATR Drugs
- FDM ATR Essential Oils
- FDM Comprehensive Organics
- FDM FTIR Polymers and Polymer Additives
- FDM FTIR Surfactants
- FDM FTIR Organic Compounds
- FDM FTIR Minerals and Inorganic Compounds
- FTIR Drugs / Canadian Forensic
- VPFTIR Organic Compounds
- FDM Forensic VPFTIR
- FDM HiRes VPFTIR



Fiveash Data Management, Inc.

608-236-9145

www.fdmspectra.com

Fax: 608-236-9170